## INSTRUCTION MANUAL

**Microscope Digital Kit**  
**Model # 44320**

# Table of Contents

Table of Contents .....	3
Introduction .....	3
Standard Accessories with your Microscope .....	4
Specifications .....	4
Magnification Table.....	5
Setting Up Your Microscope .....	5
Microscope Operation .....	6
Viewing a Specimen .....	6
Focusing & Changing Power (Magnification) .....	6
Illumination .....	6
Using the Digital Camera for Viewing and Imaging with your Microscope.....	7
Attaching your Digital Camera to your Microscope.....	11
VP-EYE Software.....	11
Care, Maintenance, and Warranty .....	13

## Introduction

Congratulations on your purchase of a Celestron microscope. Your microscope is a precision optical instrument, made of high quality materials to ensure durability and long life. It is designed to give you a lifetime of pleasure with a minimal amount of maintenance.

Before attempting to use your microscope, please read through the instructions to familiarize yourself with the functions and operations to maximize your enjoyment and usage. See the microscope diagrams to locate the parts discussed in this manual.

The microscope provides high powers from 40x to 600x. It is ideally suited for examining specimen slides of yeasts and molds, cultures, plant and animal parts, fibers, bacteria, etc. You can also examine small and thin objects at low powers such as coins, rocks, insects, various materials, etc.

With the included digital camera and the software, you can observe the magnified images or capture video or take snapshots.

The final section provides simple care and maintenance tips for you to follow to ensure that your microscope provides you with years of quality performance, usage, and enjoyment.


Figure 1

## Standard Accessories with your Microscope

• 10x—20x Zoom Eyepiece	• Digital Camera
• 4x, 15x, 30x Objective Lenses	• USB Cable
• Top & Bottom Illuminators	• Honeybee Wing, and Rock pieces, and Tweezers
• Light Diffuser	• 3 Prepared Slides --- Head of a Mosquito, Cross Section of Bamboo, Scaly Hairs of a Silverberry
• CD-ROM, Installation	• Hole (clear) Slide

## Specifications

Model # 44320	Specifications
Stage	Plain Stage with clips -- 74mm x 70mm
Zoom Eyepiece	Glass optics. Power continuous from 10x to 20x
Focuser	Coarse focus – dual knobs
Objectives	All glass optics – see magnification chart for powers
Illuminator -- Top	Pen light style. Uses 2AAA batteries (user supplied)
Illuminator -- Bottom	Uses 2AA batteries (user supplied)
Nosepiece	Triple with click stop
Camera Resolution	VGA 640x480 pixels
Weight (with batteries) and Dimensions	22oz. (624g) --- 5.25" x 3.13" x 9.75" (133mm x 79mm x 248mm)

## Magnification Table

Use the following table to determine the magnification of the different eyepiece/objective lens combination of your microscope.

Objective Lens	4x	15x	30x
10x on Zoom Eyepiece	<b>40x</b>	<b>150x</b>	<b>300x</b>
20x on Zoom Eyepiece	<b>80x</b>	<b>300x</b>	<b>600x</b>

## Setting Up Your Microscope

1. Take the Styrofoam container out of the carton.
2. Remove the tape from the Styrofoam container holding the various parts in place.
3. Carefully remove the microscope and other parts from the container and set them on a table, desk, or other flat surface.
4. Remove the plastic bag covering the microscope.
5. Remove the plastic cap from the zoom eyepiece (1).
6. Install the batteries in the top illuminator (11) which uses two AAA batteries (user supplied). See the image below (figure 2) --- Turn the top of the pen light (top illuminator) counterclockwise to take it off. Then, insert the AAA batteries into the tube with the positive ends going into the tube first. Screw the top of the pen light clockwise until tight.


**Figure 2**

7. Install the batteries for the bottom illuminator in the base (7) of the microscope which uses two AA batteries (user supplied). See the image below (3a) showing the battery compartment closed at the back of the microscope. Then, image 3b shows the batteries being installed in the compartment (see the inside of the battery compartment door to see where the positive (+) and negative (-) ends of the batteries go). Pull out on the door to open it while holding the base firmly and push firmly to close it after the batteries are installed.

**Figure 3a**

**Figure 3b**

**You are now ready to use your microscope for looking at specimen slides or small objects through the microscope zoom eyepiece! For using the digital (CMOS) camera, installation and operating instructions will follow later in this manual.**

# Microscope Operation

## Viewing a Specimen

Carefully place a specimen slide under the stage clips (9) and center the specimen directly over the hole in the center of the stage (5) – Figure 4a below shows the stage area with the hole in the center and Figure 4b shows a specimen slide centered over the hole in the stage. It will take some experimenting to place slides or objects in the center of the stage as the image you see is upside down and reversed but after some usage you will have an easy time centering. Read the sections below on Focusing, Changing Power, and Illumination before proceeding.

You are now ready to focus and view the specimen, but first you must take some precautions so as not to damage a specimen slide or valuable object. **When using the higher powers while you are focusing, make sure that the objective lens does not hit the slide or object being viewed.**

Figure 4a


Figure 4b

## Focusing & Changing Power (Magnification)

Now that the specimen slide (or object) is placed directly under the objective lens, use the focus knob (8) to focus on the specimen. Note that for very small objects, you should set them on the clear (hole) slide with a recessed hold in the center.

1. Always start with the lowest power (4x objective lens) and have the zoom eyepiece at the 10x position (all the way counterclockwise [when you are facing the microscope from the front] until it stops) so that the total power is 40x.
2. For slightly higher power, you can rotate the knurled ring on the zoom eyepiece (see **Figure 5** at the right) clockwise to obtain powers of 40x to 80x (or anywhere in between) as you continue rotating to the 20x eyepiece position. Note that you will have to refocus whenever you rotate the eyepiece to obtain a sharp focus.
3. For much higher powers, you will have to rotate the nosepiece (3) to change the objective lens to the 15x one (provides total power of 150x to 300x depending on what position you have the zoom eyepiece in) or the 30x one (provides total power of 300x to 600x). You rotate the nosepiece by holding the microscope above the nosepiece with one hand and rotate the nosepiece with the other hand until it clicks at the position. Be cautious not to let the objective lens touch the specimen slide or object when changing to higher powers – you should turn the focus knob first to lower the stage to a low position. **Note the power range of the objective lens you are using is shown on the nosepiece after it clicks into position.**
4. At the highest powers, your views will be greatly magnified but somewhat darker. The most enjoyable views can be at the lower powers which have a wider field of view and brighter illumination.

## Illumination

To get the sharpest and best views, the illumination (lighting) will have to be adjusted.

1. The top illuminator (11) is used only for solid objects (not specimen slides) so that light shines down onto the object. Push the button on the top of the pen light illuminator to turn it on/off. You can change the brightness by moving the illuminator up/down or by rotating it left or right. After some usage, you can determine the best way of adjusting the light to provide the most pleasant views.
2. The bottom illuminator (6) is used for specimen slides which shine up through the hole in the stage through the slide. You can see a close up of the illuminator in Figure 6a below. The illuminator is turned on by rotating the illuminator so the light goes up through the hole. You turn off the illuminator by rotating it down so that the mirror is on the upper part (the mirror is not useful with this microscope since you have the much better electric illuminator).
3. The light from the bottom illuminator can be increased or decreased by rotating the illuminator with very slight movements. Like with the top illuminator, you will have to experiment to provide the best lighting for the best views.

- The bottom illumination may be too bright with some specimen slides. Included with your microscope is a light diffuser which reduces the brightness and glare somewhat and can make the views sharper with a higher contrast level. In Figure 6a the diffuser is the small black piece. The diffuser fits over the bulb area by press fitting it on. Figure 6b shows what the illuminator looks like with the light diffuser in place.

Figure 6a

Figure 6b

## Using the Digital Camera for Viewing and Imaging with your Microscope

Before using your microscope for viewing on your computer screen or imaging, you first need to install the driver, software, and related files. This is a short and easy process and will only take a few minutes.

**The digital camera works with the following Microsoft Windows Operating Systems – 98/98SE/2000/ME/XP/Vista.** You also need to have a CD or DVD drive and an open USB port.


First, put the included installation CD-ROM in your computer CD/DVD drive. **Do not connect the USB cable of the imager to the computer until the driver and software are installed on your computer as the imager will not work properly.**

The following screens (left to right) will appear during the installation on most computers although variations may appear on some operating systems or even with the same operating system. The screens below are shown with an XP operating system.

The first screen to appear after you insert the CD-ROM is shown below left. You will be asked which action you wish to take to open the CD-ROM. Shown is the typical choice to make but you can make your own choice and then click “OK”.


From the next screen that appears below, double click on “VP-Eye4.0.”


Now you will double click on “Setup” using the file that is highlighted on the screen directly below.


On the screen below for the installation, please note that the Destination is “C” which the program believes is your hard drive location. If your hard drive on your computer is different, then go to the circle that says “Choose Directory” and click on it and change the “C” to the letter of your hard drive. Then, click ok.


Now click on the large circle in the center “VP EYE” and installation will begin.


During installation you will see the screen below.

The next screen you will see shows you that the software was installed properly and click “Finish”. You can check your desktop and you should see the shortcut icon “VP-EYE” and this confirms all is ok.


Go back to the second screen seen you saw originally but this time click on “Driver” as shown below.



From the screen below, click on “Setup” as shown.


Next, you will see the “InstallShield Wizard” as below and click “Next”.


Before installation of the driver is complete, you may see the screen shown below with Windows XP. If so, click on “Continue Anyway” as the program for the driver will not harm your computer.


Now you are ready to plug in the USB cable to your computer. The first screen to appear will be “Welcome to the Found New Hardware Wizard”. Put the Installation CD-ROM in to your CD/DVD drive and then click “Next” to continue.

The screen below will show the driver installation is complete. Before clicking on “Finish”, close all open files on your computer and remove the installation CD. Then, restart your computer.


The next screen indicates that the wizard is searching for the new hardware.


Before installation of the new hardware is complete, you may see the screen below with Windows XP. If so, click on “Continue Anyway” as the program for the hardware will not harm your computer.


The screen below shows that the wizard has found the hardware and is now installing the software for it.

The next screen shows that the new hardware software installation is complete. Now, remove the CD-ROM from your drive and close all open programs. Then,click “Finish”.

You should now see the screen below which prompts you to restart your computer and click “Yes”. If you do not see this screen, restart your computer anyway.

## Attaching your Digital Camera to your Microscope

After installing the driver and software, you are ready to attach your camera to the microscope.

First, you need to remove the zoom eyepiece (1) from the eyepiece tube (2) by turning it past the stop at the 10x position. Put a little pressure on the eyepiece (and one hand on the base for support) and continue turning it counterclockwise to unthread it from the eyepiece tube. Take off the protective cap from the camera. Next, thread the camera into the eyepiece tube (clockwise) until you feel it almost tight and quit. Lastly, you plug the USB cable into an open USB port on your computer.

When viewing or imaging a specimen slide or object, you can change the orientation of the image on the computer screen by rotating the camera to the position you desire – generally do this counterclockwise so you don't tighten the camera in the eyepiece tube.

Left to right below – digital camera with USB cable, eyepiece tube, camera attached to eyepiece tube, microscope with camera attached to a computer with the USB cable.

Figure 7a

Figure 7b

Figure 7c

Figure 7d

## VP-EYE Software

The software package you installed is called VP-EYE. The software allows you to observe specimen slides or objects on your computer. When you view with the camera installed, the magnification depends on the objective lens you are using and also the size of your PC monitor.

### Magnification Using Digital Camera

4x Objective using a 14" monitor ---	240x	using a 17" monitor ---	290x
15x Objective using a 14" monitor ---	850x	using a 17" monitor ---	1030x
30x Objective using a 14" monitor ---	1550x	using a 17" monitor ---	1880x

The VP-EYE software contains many applications capturing and processing your photos (snapshots) or videos.


All applications can be started from the Application Panel. It can be started by clicking on the VP-EYE icon or by accessing the programs on your computer.

**Note: The software package may not have some of the programs listed due to the version of this software package.**


From the Application Panel (Settings) you can change the video drivers, language, start up menu, auto-run features, etc. On the Application Panel are various buttons (see screen below) and just put the cursor over them to see what they are. Each screen has a help button to give you detailed information about each. The main function buttons (Photo Bank, Photo Function, Video Bank, and Video Function) are shown below. You can print out a short instruction sheet from the installed software.

If you need additional information about the software, go to their website at <http://www.mmedia.com.tw>

Application Panel to the right ----

**Photo Bank**

**Photo Function**

**Video Bank**


**Video Function**

Below are a few snapshot images taken by a young teenager on his first attempts after installing the driver and software.


**Winter Jasmine Leaf w/4x Objective**


**Winter Jasmine Leaf w/15x Objective**


**U.S. Penny w/4x Objective**

**Year Tilia Stem w/4x Objective**

**Year Tilia Stem w/15x Objective**

**Rock w/4x Objective**

## Care, Maintenance, and Warranty

Your Celestron microscope and digital camera are precision optical instruments and should be treated with care at all times. Follow these care and maintenance suggestions and your microscope will need very little maintenance throughout its lifetime.

- When you are done using your microscope, remove any specimens left on the stage.
- Turn off the top and bottom illuminators when you are done using the microscope.
- If you will not be using your microscope for a long period of time, remove the batteries from the top and bottom illuminators.
- Always place the dust cap over the eyepiece and the camera sensor when not in use or when being stored.
- Store the microscope in a dry and clean place.
- Be very careful if using your microscope in direct sun light to prevent damage to the microscope or your eyes.
- Never point the sensor on the camera towards the sun or the camera can be damaged and cease working.
- When moving your telescope, carry it by the “arm” with one hand.
- Clean the outside surfaces with a moist cloth.
- Never clean optical surfaces with cloth or paper towels as they can scratch optical surfaces easily.
- Blow off dust with a camels hair brush or an air blower from optical surfaces.
- To clean fingerprints off of optical surfaces, use a lens cleaning agent and lens tissue available at most photo outlets and when cleaning do not rub in circles as this may cause streaks or scratches to occur.
- Never disassemble or clean internal optical surfaces. This should be done by qualified technicians at the factory or other authorized repair facilities.
- When handling glass specimen slides, use care as the edges can be sharp.

### Replacing the bottom illuminator bulb

Your microscope comes with one spare bulb for the bottom illuminator. See Figure 6a and you will see the white portion of the illuminator. Turn the small cap on the white housing counterclockwise while pushing in to remove it. Push in and rotate counterclockwise the bulb to be replaced. Take the spare bulb and do the reverse of the above steps.

### Warranty

Your microscope has a two year limited warranty. Please see the Celestron website for detailed information on all Celestron microscopes at [www.celestron.com](http://www.celestron.com).

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Torrance, California 90503 U.S.A.